

MORTUARY TROLLEY

- Mild tubular steel frame – epoxy coated.
- Stainless steel top section.
- With push handles at each end of trolley.
- Unit is supplied complete with removable canopy frame with canvas cover.

Castors: 2 x 600mm Fixed directional wheels mounted at centre of trolley.
2 x 150mm Swivel Castors

Dimensions:

Length: 2180mm
Width: 600mm
Height: 730mm
Weight: 49 kg
Volume: 0.95 m³

A4-087-ES6



MORTUARY TROLLEY

- Mild steel construction – epoxy coated.
- Stainless steel or epoxy top.
- **Option:** Canopy frame with canvas.

Castors: 4 x 150mm swivel castors (2 with brakes).

Dimensions:


Length: 2180mm
Width: 600mm
Height: 730mm
Weight: 30 kg
Volume: 0.95 m³

M103

S/STEEL TOP AND M/STEEL HANDLES

M103A

M/STEEL



MORTUARY STACKER TROLLEY

HYDMORT

ELECTRIC MORTUARY TROLLEY

ELECMORT

ELECTRIC MORTUARY TROLLEY

- Lifts from a height of 300mm to 1800mm
- Easily moved on castors, rear castors fitted with brake
- Complete s/steel construction
- Hydraulic or Electric lifting action
- Max weight capacity 180kg
- Double action manual hydraulic pump for quick lifting
- Controlled release for easy and save lowering
- Accepts all standard 500mm wide mortuary body trays
- Body tray L2265 x W500mm

Dimensions:

Length: 2000 mm
Width: 650mm
Height: 300 to 1800mm
Weight: 90kg

Electric / Hydraulic model

- 12 V DC Electric pump.
- 220v Battery charger
- 100Ah Battery on board
- High speed lift and lower.
- Overnight charge (8 hours) allows 20 full height lifts.
- Easy one hand / one lever operation



**SAVE AND EASY WAY TO LOAD
CADAVERS INTO MORTUARY CABINETS
OR COLD ROOM RACKS.**

ECONOMY / STANDARD / MEDICAL 3, 4 OR 6 TIER

- Chromadek interior and exterior
- S/steel doors and front
- S/steel inner floor tray and drain
- S/steel body trays
- Galvanised lower frame and levelling feet
- Nylon rollers on s/steel pins
- Rear mounted removable refrigerator system
- Add 380mm to the length for rear mounted refrigeration

Dimensions:	3 Tier	4 Tier	6 Tier
Length:	2900mm	2900mm	2900mm
Width:	890mm	1780mm	1780mm
Height:	1900mm	1465mm	1900mm
Weight:	1000kg	1500kg	2000kg



COLD ROOM RACKS

- Racks and body trays can be installed depending on room size.



BODY TRAY

- S/steel body tray available in standard size of 550mm or Obese 650mm wide.
- S/steel forensic body tray available in standard size of 550mm or Obese 650mm wide.
- Forensic body tray has an outlet for fluids

BTRAY550

BTRAY650

FBTRAY550

FBTRAY650

