

MEDCOR MEDICAL EQUIPMENT (PTY) LTD t/a



M E D I T E K *hemco*
MEDICAL EQUIPMENT SALES

PRODUCT CATALOGUE

ACCESSORIES

- > STEPS
- > DRIP STANDS
- > ARM STANDS
- > KICK-ABOUT
- > BEDSCREENS
- > RACKS
- > BABY BASSINET CRIB
- > X-RAY VIEWING BOX
- > HOLDERS
- > COLLAPSIBLE CRADLE



DOUBLE BED STEP

M36

DESCRIPTION

- Tubular steel frame – epoxy finish
- Two steps with rubber matting and aluminium edging
- Four feet with rubber protectors

DIMENSIONS

Length	500mm
Width	490mm
Height	325mm



DOUBLE BED STEP - S/STEEL

M36SS

DESCRIPTION

- Stainless steel frame
- With aluminium chequered plate tops
- Four rubber feet

DIMENSIONS

Length	500mm
Width	490mm
Height	325mm



SINGLE BED STEP

M35

DESCRIPTION

- Tubular steel construction with epoxy coated finish
- ABS rubber top with aluminium edging
- Feet mounted on rubber protectors

DIMENSIONS

Length	530mm
Width	335mm
Height	160mm



SINGLE BED STEP S/STEEL

M35SS

DESCRIPTION

- Tubular S/steel frame
- With aluminium chequered plate top
- Feet mounted on rubber protectors

DIMENSIONS

Length	530mm
Width	335mm
Height	160mm



DRIP STAND

M18A

DESCRIPTION

- Tubular steel frame epoxy coated
- Five star nylon base on castors
- Height adjustable
- With two s/steel hooks for drip bags

DIMENSIONS

Diameter	450mm (base)
Height	1300-2100mm



DRIP STAND S/STEEL

M18ASS

DESCRIPTION

- Tubular steel frame Stainless steel
- Five star s/steel base on castors
- Height adjustable
- With 2 hooks for drip bags

DIMENSIONS

Diameter	450mm (base)
Height	1300-2100mm



DRIP STAND

M18AE

DESCRIPTION

- Tubular steel frame epoxy coated
- Five star epoxy coated steel base on castors
- Height adjustable
- With s/steel 2 hooks T rod

DIMENSIONS

Diameter 450mm (base)
Height 1300-200mm



DRIP STAND WITH HOLDER

M18AHOLD

DESCRIPTION

- Tubular steel frame epoxy coated
- Five star epoxy coated base on castors
- Height adjustable
- With 2 hooks for drip bags
- With holder

DIMENSIONS

Diameter 450mm (base)
Height 1300-2100mm



DRIP STAND

M17

DESCRIPTION

- Tubular steel frame work epoxy coated finish
- Three star base on castors
- Height adjustable
- Stainless steel hooks for drip bags

DIMENSIONS

Diameter 450mm (base)
Height 1300-2100mm



CLAMP ON DRIP ROD

M16

DESCRIPTION

- Round bar mild steel epoxy coated construction
- Two stainless steel hooks for IV bags
- Rod length 1500mm – fully adjustable
- Clamp on facility for square tubular frame

DIMENSIONS

N/A



BLOOD DONOR ARM STAND WITH HOLDER

BLOODPATH

DESCRIPTION

- Round tubular steel frame with epoxy coated finish
- Upholstered armrest
- S/steel holder for kidney dish size 250mm
- Adjustable Height
- 3 Star base

DIMENSIONS

Length	350mm
Width	350mm
Height	690mm



BLOOD DONOR ARM STAND

BLOOD

DESCRIPTION

- Round tubular steel frame with epoxy coated finish
- Upholstered armrest
- Adjustable Height
- 3 Star base

DIMENSIONS

Length	350mm
Width	350mm
Height	690mm



KICK-ABOUT

M53A

DESCRIPTION

- Steel frame - epoxy coated finish
- Optional: Stainless steel 12L bucket
- **Optional:** Lid for bucket

DIMENSIONS

Diameter	385mm
Height	330mm



KICK-ABOUT

M53

DESCRIPTION

- Round tubular steel frame with epoxy coated finish
- **Optional:** Stainless steel bowl

DIMENSIONS

Diameter	385mm
Height	330mm



BEDSCREEN - MOBILE

M29	(3 Section)
M29A	(2 Section)
M29B	(4 Section)

DESCRIPTION

- Tubular steel frame with epoxy coated finish
- Folding screens with non-transparent washable curtains
- Mobile unit

DIMENSIONS

Panel	600mm (Unit Width)
Height	1500mm (Panel) / 1800mm (Overall)



BEDPAN RACK

M27

DESCRIPTION

- Grade 304 Stainless steel construction
- Capacity - 5 bedpans
- With removable grade 304 s/steel drip tray
- Wall mounted unit

DIMENSIONS

N/A



URINAL RACK

M28

DESCRIPTION

- Grade 304 Stainless steel construction
- Capacity - 5 urinals
- With removable grade 304 s/steel drip tray
- Wall mounted unit

DIMENSIONS

URINAL RACK



BABY BASSINET CRIB

CRIBO – Crib with water outlet
 CRIB – Crib without water outlet

DESCRIPTION

- Clear Ampalite construction
- With or without water outlet
- **Optional:** Mattress with fluid proof cover

DIMENSIONS

Length	740mm
Width	440mm
Height	260mm





X-RAY VIEWING BOX - LED

- MEC49 - Small
- MEC49 - Medium
- MEC49 - Large

DESCRIPTION

- Aluminium frame
- Uniform lighting
- Perspex screen and slide

DIMENSIONS

(VISUAL DIMENSIONS)

	Small	Medium	Large
Length	370mm	730mm	1080mm
Height	440mm	440mm	440mm(Overall)
Depth	110mm	110mm	110mm



CHART AND X-RAY HOLDER

MEC46

DESCRIPTION

- Sheet metal construction with epoxy coated finish
- With 2 compartments to hold documentation
- Facility to hook over bed end

DIMENSIONS

Height	350mm
Width	385mm
Depth	40mm



HOLDER - URINAL BAG

M85

DESCRIPTION

- Wire construction – epoxy coated or stainless steel
- Hook on facility to mattress section of bed

DIMENSIONS

Length	500mm
Width	200mm



COLLAPSIBLE LEG & BODY CRADLE

M19 - Large

M20 - Medium

M21 - Small

DESCRIPTION

- Mild steel rod with epoxy coated finish
- Fold down for easy storage


DIMENSIONS

	Small	Medium	Large
Width -	450mm	660mm	660mm
Length -	535mm	610mm	800mm
Height -	300mm	435mm	450mm



REQUEST A QUOTE

 carla@meditek.co.za

 021 553 1044 | 084 422 6603

 www.meditek.co.za