

MEDCOR MEDICAL EQUIPMENT (PTY) LTD t/a



M E D I T E K *hemco*
MEDICAL EQUIPMENT SALES

PRODUCT CATALOGUE

BEDS & COTS

- > BEDS
- > COTS
- > BASSINET
- > NURSE DIVAN



BED FULLY ELECTRIC

M2ELECPC

DESCRIPTION

- Mild steel construction with epoxy coated finish
- Electric operated Height, Backrest and Knee break sections
- Hand control unit
- Including guardrails & IV Pole
- Max user weight 250kg
- **Supplied in knock down form - Assembly required**

DIMENSIONS

Length	1920mm
Width	900mm
Height	400 – 730mm



WARD BED ADJUSTABLE HEIGHT

M2PCH

DESCRIPTION

- Steel construction epoxy coated
- Mild steel mattress seat
- Manual (winding) operated Height, Backrest and Knee break sections
- User weight capacity 250kg
- Included: Guardrails & IV Pole
- **Supplied in knock down form - Assembly required**

DIMENSIONS

Length	2000mm (Mattress surface)
Width	900 mm
Height	360 – 710mm



BED FIXED HEIGHT

M2PC

DESCRIPTION

- Steel construction epoxy coated
- Mild steel mattress seat
- Fixed height
- Winding operated Backrest and Knee break sections
- User weight capacity 250kg
- Including: Guardrails & IV Pole
- **Supplied in knock down form - Assembly required**

DIMENSIONS

Length	2000mm (Mattress surface)
Width	900 mm
Height	500mm



NURSE DIVAN

M5 – SINGLE SECTION

M5A – GASLIFT

M5B - WINDING

DESCRIPTION

- Mild steel mattress section – epoxy coated
- Backrest Options: Flatbed / Gaslifted / Winding
- Removable head and foot ends

DIMENSIONS

Length	2000mm
Width	915mm
Height	530mm (mattress base)



NURSES DIVAN - STATIONARY

MEC39

DESCRIPTION

- Single section bed with weld mesh mattress section
- Mild steel construction – epoxy coated finish
- Mounted on 4 legs with rubber protectors
- Supplied in knock-down form

DIMENSIONS

Length	2038mm
Width	915mm
Height	530mm



CHILDREN'S COT – MOBILE SINGLE SECTION

A1-007-002

DESCRIPTION

- Round tube steel construction with epoxy coated finish
- Wire mesh mattress seat epoxy coated.
- Cot sides of drop-down type
- Tubular head and foot end of equal height
- **Optional:** Mattress

DIMENSIONS

Length	1370mm
Width	762mm
Height	685mm (To mattress seat)



COT - ADULT

- M10 - Single section
- M10C - With Gas lift Backrest

DESCRIPTION

- Tubular steel construction with weld mesh mattress seat – epoxy coated
- Two sliding sides with automatic locking device
- **Optional:** Mattress

DIMENSIONS

Length	2000mm
Width	915mm
Height	685mm (To mattress seat) 1200mm (To top of headboard)



ARMINIA ERGO ELECTRIC 4 SECTION

ARMINIA

DESCRIPTION

- Steel construction
- Slatted mattress seat
- Electric operated Height, Backrest and Knee break sections
- Patient hand control unit
- Beech wooden panelled head and foot boards
- Collapsible safety sides in beech wooden finish
- Lifting pole
- User weight capacity 145kg
- **Note: Assembly required**

DIMENSIONS

Length	2240mm (Outer) / 2080mm (Inner)
Width	1020mm (Outer) / 900mm (Inner)
Height	400 - 800mm
Backrest angle	70 degrees
Thigh angle	24 degrees



MOBILE BASSINET

M9A

DESCRIPTION

- Tubular steel framework with sheet metal blanket tray – epoxy coated
- With Trendelenberg position
- With plastic crib no outlet
- **OPTIONAL:** Crib with outlet
- **OPTIONAL:** With Cordoba covered mattress


DIMENSIONS

Length	750mm
Width	355mm
Height	800mm



REQUEST A QUOTE

 carla@meditek.co.za

 021 553 1044 | 084 422 6603

 www.meditek.co.za