

MEDCOR MEDICAL EQUIPMENT (PTY) LTD t/a



M E D I T E K *hemco*  
MEDICAL EQUIPMENT SALES

# PRODUCT CATALOGUE

## CARTS & CABINETS

---

- > DRUG CABINETS
- > INSTRUMENT CABINETS
  - > CARTS
- > BED END TROLLEYS
  - > TROLLEYS
- > MOBILE STERIPAK



## DRUG CABINET - SCHEDULE 7

A0-T45B-E00

### DESCRIPTION

- Sheet metal construction - epoxy coated finish
- Cabinet with sloping top
- Special safety locks and keys
- Double lockable doors with internal lockable schedule 7 compartment
- Wall mounted unit
- **Optional:** Poison / Schedule 7 sticker

### DIMENSIONS

Height	711mm
Width	7110mm
Depth	280mm



## DRUG CABINET - SCHEDULE 5 & 6

M31 (Single Door - Small)  
M31A (Single Door - Large)

### DESCRIPTION

- Sheet metal construction with epoxy coated finish
- Lockable doors with special safety lock and key
- With one internal shelf
- Wall mounted
- **Optional:** Schedule 5 & 6 sticker

### DIMENSIONS

<b>Small</b>	
Height	460mm
Width	360mm
Depth	230mm
<b>Large</b>	
Height	610mm
Width	460mm
Depth	305mm



## DRUG CABINET - SCHEDULE 5 & 6

M32 (Double Door - Small)  
M32A (Double Door - Large)

### DESCRIPTION

- Sheet metal construction with epoxy coated finish
- Lockable double door compartments with special safety lock and key each
- With one internal shelf
- Wall mounted
- **Optional:** Schedule 5 & 6 sticker

### DIMENSIONS

<b>Small</b>	
Height	460mm
Width	360mm
Depth	230mm
<b>Large</b>	
Height	610mm
Width	460mm
Depth	305mm



## INSTRUMENT CABINET - DBL DOOR

M80

### DESCRIPTION

- Sheet metal construction with epoxy coated finish
- Mounted on adjustable feet
- Two lockable glass doors side by side cover the full length of the cabinet
- The cabinet is supplied with 6 adjustable glass shelves
- **Optional:** Unit can be supplied on castors

### DIMENSIONS

Length	970mm
Depth	500mm
Height	1970mm



## INSTRUMENT CABINET

M81A

### DESCRIPTION

- Sheet metal construction with epoxy coated finish
- Mounted on adjustable feet
- Upper section with 2 lockable sliding doors
- 2 Adjustable glass shelves
- Lower section with 2 lockable metal doors
- With internal sheet metal removable shelf
- **Optional:** Unit can be supplied on castors

### DIMENSIONS

n/a



## EMERGENCY CART

MX34EMG

### DESCRIPTION

- Sheet metal construction with non toxic epoxy coated finish in Red
- ABS top
- 5 drawers (2 x 3" / 2 x 6" / 1 x 9")
- Central locking system
- 3" Drawer divider set
- 6" Drawer divider set
- Slide out working surface
- Push handle situated on one side
- IV Pole
- VDS Defibrillator shelf
- Oxygen cylinder holder & Side rail
- CPR Board
- Moulded medication tray
- Sharps bin bracket for rectangular or round bins
- **Optional extra:** Accessories are available on request
- **Note:** MX37EMG – Emergency cart 6 drawer also available

### DIMENSIONS

Length	690mm
Depth	518mm
Height	1006mm



## PROCEDURE CART

MX34PRO

### DESCRIPTION

- Sheet metal construction with non toxic epoxy coated finish in White
- ABS top
- 5 drawers (2 x 3" / 2 x 6" / 1 x 9")
- Central locking system
- Slide out working surface
- 3" Drawer divider set
- 6" Drawer divider set
- Push handle situated on one side
- Sharps bin bracket for rectangular or round bins
- **Optional:** Accessories are available on request

### DIMENSIONS

Length	690mm
Depth	518mm
Height	1006mm



## ANESTHESIA CART

MX37ANE

### DESCRIPTION

- Sheet metal construction with non toxic epoxy coated finish in Blue
- 6 Drawers (3 x 3" / 2 x 6" / 1 x 9")
- Central locking system
- Top drawer for Narcotics - lockable
- Bottom drawer incl. dividers
- Slide out working surface
- Push handle situated on one side
- Sharps bin bracket for rectangular or round bins
- **Optional:** Accessories are available on request

### DIMENSIONS

Height	1085mm
Width	690mm
Depth	518mm



## PAEDIATRIC CART

MXPAED

### DESCRIPTION

- Sheet metal construction with non toxic epoxy coated finish in pastel Light blue
- Child friendly Sea Sticker design
- ABS top
- 5 drawers (2\*3" / 2\*6"/1\*9")
- Central locking system
- Slide-out work surface
- Push handle situated on one side
- 3" Fixed divider Set
- 6" Fixed divider set
- Sharps bin bracket for rectangular or round bins
- **Optional:** Available on request

### DIMENSIONS

Height	1006mm
Width	690mm
Depth	518mm





## PAEDIATRIC EMERGENCY CART

**BROSELOW**

### DESCRIPTION

- Sheet metal construction with non toxic epoxy coated finish
- ABS top
- 9 Colour coded drawers (8\*3"/1\*6") for emergency Paediatric care
- Central locking system
- Break-away lock for individual drawers
- Slide-out work surface
- Push handle situated on one side
- **Optional:** Sharps bin-bracket only (Round/ Rectangular) including side-bar mounting rail VPR-2
- **Optional:** available on request (Does not include any drawer divider sets)

### DIMENSIONS

Height	1085mm
Width	690mm
Depth	518mm



## MEDICATION CART - DOUBLE SIDED

**FC3441DXL**

### DESCRIPTION

- Sheet metal construction with non toxic epoxy coated finish
- ABS top
- 4 Medication Drawers (4 x 3" (20 x 3" bins) / 4 x 6" (20 x 6" bins) / 1 x 9")
- Central locking system
- Push handles situated on both sides
- **Optional:** Accessories are available on request

### DIMENSIONS

Height	1090mm
Width	790mm
Depth	490mm



## ICU BED END TROLLEY

**BEDEND105**

### DESCRIPTION

- Mild steel epoxy coated construction
- 5 Drawers
- 1 Cupboard section with internal shelf
- Formica wooden writing service top with aluminium edging at the bottom
- Shelf situated underneath the writing surface
- With holder for pens etc. situated at top right corner
- X Ray chart holder on the side of the unit
- Edging around 3 sides of shelf under Formica top

### DIMENSIONS

Length	1150mm
Width	600mm
Height	1000mm



## BED END TROLLEY

BEDEND104

### DESCRIPTION

- Mild steel epoxy coated construction
- 2 Drawers
- Formica wooden writing service top with aluminium edging at the bottom
- Shelf situated underneath the writing surface
- X Ray chart holder on the side of the unit

### DIMENSIONS

N/A



## NEONATAL BED END TROLLEY

BEDEND106

### DESCRIPTION

- Mild steel epoxy coated construction
- 1 Drawer
- Formica wooden writing service top with aluminium edging at the bottom
- Shelf situated underneath the writing surface
- 1 Pull out wire basket

### DIMENSIONS

Length	600mm
Width	600mm
Height	1100mm (front), 1250mm (back)



## ANAESTHETIC TROLLEY

A2-026C-ES5

### DESCRIPTION

- Mild steel construction with epoxy coated finish
- Grade 18/12 stainless steel top
- Cabinet has 6 drawers and 3 shelves
- **Optional:** Top drawer fitted with a lock

### DIMENSIONS

Length	1070mm
Width	460mm
Height	1100mm



## ANAESTHETIC TROLLEY – DRAWERS ONLY

A2-026CAB

### DESCRIPTION

- Mild steel construction with epoxy coated finish & s/steel top
- Grade 18/12 stainless steel top
- Cabinet has 6 drawers

### DIMENSIONS

Length	690mm
Width	460mm
Height	1100mm



## EMERGENCY CART

M79B

### DESCRIPTION

- Mild steel epoxy coated construction
- S/steel top with guardrail around 3 sides
- Five lockable drawers
- Complete with oxygen cylinder holder, drip rod, writing shelf and CPR board

### DIMENSIONS

Length	900mm
Width	450mm
Height	850mm



## CABINET – MOBILE STERIPAK

A4-077A-EE4

### DESCRIPTION

- Mild steel epoxy coated construction
- Lockable doors
- Two adjustable shelves (internal)
- Push handle mounted on one side
- **Optional:** S/steel shelf below cupboard section, push handle on opposite side of cabinet


### DIMENSIONS

Length	870mm
Width	410mm
Height	1200mm



## REQUEST A QUOTE

 [carla@meditek.co.za](mailto:carla@meditek.co.za)

 021 553 1044 | 084 422 6603

 [www.meditek.co.za](http://www.meditek.co.za)