

MEDCOR MEDICAL EQUIPMENT (PTY) LTD t/a



M E D I T E K *hemco*
MEDICAL EQUIPMENT SALES

PRODUCT CATALOGUE

GENERAL TROLLEYS

- > INSTRUMENT TROLLEYS
 - > DRESSING TROLLEYS
 - > ANAESTHETIC TROLLEY
 - > BOWL STAND
 - > MEDICINE TROLLEY
 - > POP TROLLEY
 - > TEA/FOOD TROLLEYS
 - > MOBILE UTILITY TROLLEY
 - > OXYGEN CYLINDER STAND & TROLLEYS
 - > SWAB RACK TROLLEY
 - > MOBILE BEDPAN AND URINAL STAND



INSTRUMENT TROLLEY

- M55 – Large
- M56 – Medium
- M57 - Small

DESCRIPTION

- Tubular steel frame – epoxy coated finish
- Two stainless steel shelves - removable

DIMENSIONS

	L	M	S
Length	915mm	600mm	450mm
Width	450mm	450mm	450mm
Height	860mm	860mm	860mm



INSTRUMENT TROLLEY S/STEEL

- M55SS – Large
- M56SS – Medium
- M57SS - Small

DESCRIPTION

- Complete S/steel construction
- Two stainless steel shelves - removable

DIMENSIONS

	L	M	S
Length	915mm	600mm	450mm
Width	450mm	450mm	450mm
Height	860mm	860mm	860mm



DRESSING TROLLEY

- M51 - Medium
- M52 - Small
- M52b - Large

DESCRIPTION

- Tubular steel frame with epoxy coated finish
- Two stainless steel shelves
- Guard rails on top and bottom shelves

DIMENSIONS

	M	S	L
Length	610mm	455mm	915mm
Width	455mm	455mm	455mm
Height	860mm	860mm	860mm



DRESSING TROLLEY S/STEEL

- M51SS - Medium
- M52SS - Small
- M52bSS - Large

DESCRIPTION

- Complete S/steel construction
- Two stainless steel shelves
- Guard rails on top and bottom shelves

DIMENSIONS

	M	S	L
Length	610mm	455mm	915mm
Width	455mm	455mm	455mm
Height	860mm	860mm	860mm



ANAESTHETIC TROLLEY SMALL

M59

DESCRIPTION

- Tubular steel frame with epoxy coated finish
- S/steel top with guard rail around 3 sides
- Single drawer (without lock)
- With stainless steel bottom shelf
- **Optional:** With lockable drawer
- **Optional:** Complete s/steel unit

DIMENSIONS

Length	455mm
Width	455mm
Height	860mm



ANAESTHETIC TROLLEY LARGE

M50

DESCRIPTION

- Tubular steel frame with epoxy coated finish
- S/steel top with guard rail around 3 sides
- Double drawer (without lock)
- With stainless steel bottom shelf
- **Optional:** Lockable drawers
- **Optional:** Complete s/steel unit

DIMENSIONS

Length	915mm
Width	458mm
Height	850mm



SINGLE BOWL STAND

M44

DESCRIPTION

- Tubular mild steel construction epoxy coated
- **Optional:** 1x Bowl 335mm

DIMENSIONS

Diameter	355mm
Height	900mm



SINGLE BOWL STAND S/STEEL

M45

DESCRIPTION

- Complete tubular S/steel construction
- **Optional:** 1x Bowl 335mm

DIMENSIONS

Diameter	355mm
Height	900mm



DOUBLE BOWL STAND

M46A

DESCRIPTION

- Square and round tubular mild steel frame with epoxy coated finish
- **Optional:** 2x 350mm St/Steel Bowls

DIMENSIONS

Length	800mm
Width	355mm
Height	840mm



DOUBLE BOWL STAND S/STEEL

M46

DESCRIPTION

- Square and round tubular s/steel frame
- **Optional** – 2x 350mm S/Steel Bowls

DIMENSIONS

Length	800mm
Width	355mm
Height	840mm



MOBILE MEDICINE TROLLEY

M49

DESCRIPTION

- Round tubular steel sub-frame with epoxy coated finish
- With stainless steel bottom shelf
- Cabinet – Sheet metal construction with epoxy finish
- With internal fold up writing shelf
- With sloping front
- Push handles allocated on both sides of unit
- Lockable double doors
- **Optional:** Holder for stainless steel bowl and dish
- **Optional:** S/steel round bowl and square s/steel dish
- **Optional:** Pull out writing shelf allocated below cabinet

DIMENSIONS

	CABINET	OVERALL
Length	600mm	600mm
Width	470mm	470mm
Height	500mm	1250mm



POP TROLLEY

M83 (Mild steel frame / SS Shelves)
M83A (Complete S/steel trolley)

DESCRIPTION

- Tubular steel construction with epoxy coated finish
- With two stainless steel shelves, top shelf with guardrail around 3 sides
- With holder for bucket & bowl
- **Optional:** Complete s/steel trolley (M83A)
- **Optional:** Bucket, Bowl

DIMENSIONS

Length	915mm
Width	450mm
Height	860mm



TEA/FOOD TROLLEY - 3 TIER

A2-063-EE3

DESCRIPTION

- Tubular steel frame epoxy coated
- Three epoxy coated shelves

DIMENSIONS

Length	914mm
Width	437mm
Height	1054mm (75MM CASTORS)



TEA/FOOD TROLLEY S/STEEL - 3 TIER

A2-063-SS3

DESCRIPTION

- Complete tubular S/steel construction
- Three s/steel shelves

DIMENSIONS

Length	914mm
Width	437mm
Height	1054mm (75MM CASTORS)



TEA/FOOD TROLLEY

M40 2-TIER

M41 3-TIER

DESCRIPTION

- Square mild steel frame – epoxy coated finish
- Two stainless steel shelves
- Push handles both sides
- **Optional:** Overall stainless steel

DIMENSIONS

Length	1100mm
Width	460mm
Height	810mm



MOBILE UTILITY TROLLEY S/STEEL

M88 (Small / Medium / Large)

DESCRIPTION

- Stainless steel construction
- With 2 stainless steel shelves
- With one push handle
- **Optional:** Additional extra stainless steel shelf

DIMENSIONS

	S	M	L
Length	600mm	700mm	900mm
Width	400mm	500mm	600mm
Height	800mm	800mm	800mm



MEDICINE TROLLEY S/STEEL

M89M - Medium

M89L - Large

DESCRIPTION

- Stainless steel construction
- 2 stainless steel shelves
- With 1 drawer
- Holder for bowl
- One push handle
- **Optional:** Stainless steel bowl
- **Optional:** 2 Drawers

DIMENSIONS

	M	L
Length	700mm	900mm
Width	500mm	600mm
Height	800mm	800mm



MEDICINE TROLLEY

MEC7 – Epoxy

MEC7A – S/Steel

DESCRIPTION

- Tubular steel construction – epoxy or stainless steel
- Two shelves s/steel
- Guardrails around all four sides of shelves
- **Optional:** Bottle holder

DIMENSIONS

Length	750mm
Width	400mm
Height	860mm



OXYGEN CYLINDER STAND

M30

DESCRIPTION

- Mild steel construction with epoxy coated finish

DIMENSIONS

Diameter	400mm (Overall)
Height	780mm



OXYGEN CYLINDER TROLLEY

M33

DESCRIPTION

- Mild tubular steel construction with epoxy coated finish
- Trolley can be used for various lengths of cylinders

DIMENSIONS

Length	480mm
Width	355mm
Height	940mm



OXYGEN CYLINDER STAND

A2-070-E03

DESCRIPTION

- Mild steel tubular construction with epoxy coated finish
- Cylinder can be clamped by means of vice screws

DIMENSIONS

Diameter	540mm (Overall)
Height	670mm



SWAB RACK TROLLEY

M34

DESCRIPTION

- Complete tubular stainless steel construction
- Equipped with 3 racks
- Drip tray s/steel
- **Optional:** Note board

DIMENSIONS

Length	1580mm
Width	500mm
Height	1120mm



MOBILE BEDPAN AND URINAL STAND

A2-038-SS3

DESCRIPTION

- Stainless steel construction
- Facility to carry 6 bedpans and 12 urinals
- 2 Stainless steel drip trays, one under the bedpans and one under the urinals

DIMENSIONS

Length	1077mm
Width	440mm
Height	1255mm



PHARMACY MEDICINE TROLLEY

M90

DESCRIPTION

- Tubular steel frame with melamine shelf
- Four lockable melamine drawers
- Chrome push handles both sides
- Bottom melamine shelf


DIMENSIONS

Length	800mm
Width	400mm
Height	1150mm



REQUEST A QUOTE

 carla@meditek.co.za

 021 553 1044 | 084 422 6603

 www.meditek.co.za