

MEDCOR MEDICAL EQUIPMENT (PTY) LTD t/a



M E D I T E K *hemco*
MEDICAL EQUIPMENT SALES

PRODUCT CATALOGUE

MATTRESSES

- > STANDARD RIPPLE
FOAM
- > ANTI BEDSORE
HYPERFOAM
- > EGGBOX MATTRESS
OVERLAY



STANDARD RIPPLE FOAM MATTRESS


PCRIPMAT

DESCRIPTION

- High density foam
- Fluid proof 2 way stretch cover
- Anti flammable, anti bacterial
- Air permeable
- Pressure relief
- Zip covered with flap
- Max user weight 110kg

DIMENSIONS

Length	2000mm
Width	900mm
Thickness	140mm



ANTI BEDSORE HYPERFOAM 2 CLINIC


MATH2C

DESCRIPTION

- Up to stage 4 anti-bedsore foam mattress
- Applicable for pain therapy
- Fluid proof breathable 4 way stretch cover
- Firm edges provide stability
- Zip covered with flap
- Max user weight 150kg

DIMENSIONS

Length	2000mm
Width	900mm
Thickness	150mm



ANTI BEDSORE HYPERFOAM CC

MATHYPCC

DESCRIPTION

- Up to stage 1 anti bedsore foam mattress
- Applicable for pain therapy
- Fluid proof breathable cover—4 way stretch
- Firm edges provide stability
- X ray translucent
- Zip covered with flap
- Max user weight 140kg

DIMENSIONS

Length	2000mm
Width	900mm
Thickness	140mm



EGGBOX MATTRESS OVERLAY

EGGBOS – Single bed
EGGBOD – Double bed
EGGBDQ – Queen bed
EGGBOK – King size bed


DESCRIPTION

- Helps prevent and relieve bedsores
- High density foam with egg box shape
- Available in: Single, 3/4, Double, Queen and King sizes (2×Single)



REQUEST A QUOTE

 carla@meditek.co.za

 021 553 1044 | 084 422 6603

 www.meditek.co.za