

MEDCOR MEDICAL EQUIPMENT (PTY) LTD t/a



M E D I T E K *hemco*  
MEDICAL EQUIPMENT SALES

# PRODUCT CATALOGUE

## PATIENT TRANSFER

---

- > TROLLEYS
- > STRETCHERS
- > SPINE BOARD
- > HEAD IMMOBILIZER
- > SPIDER STRAP  
ADJUSTABLE



## PATIENT TROLLEY

A2-086-ES5

### DESCRIPTION

- Tubular steel construction with epoxy coating
- Mattress seat section stainless steel
- Chrome plated push handles at each end of trolley
- 2 Drip rod fittings
- With drip rod
- With blanket tray
- With oxygen cylinder cage
- With restraining straps along length of trolley
- **Optional:** Guardrails
- **Optional:** With adjustable ratchet backrest

### DIMENSIONS

Length	2070mm
Width	620mm
Height	780mm



## RECOVERY TROLLEY

A1-016-EE6

### DESCRIPTION

- Epoxy powder coated frame
- Forward and reverse Trendelenburg operated by means of a handle
- The trolley has full length drop down sides and removable head & foot ends
- Provision for a drip rod is included at each corner of the trolley
- A drip rod is included along with a drip rod holder
- An oxygen cylinder holder and blanket tray are secured to the base of the unit
- **Optional:** Mattress

### DIMENSIONS

Length	2050mm
Width	700mm
Height	750mm



## PAEDIATRIC TRAUMA TROLLEY

AC-116-PAED

### DESCRIPTION

- Tubular construction epoxy coated
- Removable head and foot ends
- X-Ray Facility **ON REQUEST**
- With Trendelenburg position
- With Oxygen cylinder holder
- With 4x driprod clamps (2x on each side) and IV pole
- With tray for accessories
- Including fluidproof thick PVC mattress
- Including fold-down guardrails
- Maximum weight capacity: 100kg

### DIMENSIONS

Length	1630mm
Width	680mm



## AMBULANCE STRETCHER


DDC-3

### DESCRIPTION

- Aluminium construction
- Mattress surface is a 50mm thick soft foamed cushion for patient comfort
- Foldaway mechanical legs by means of handles on both sides
- Height adjustable
- Max weight capacity 159kg

### DIMENSIONS

Length	1930mm
Width	570mm



## SPINE BOARD

DDJ-6B

### DESCRIPTION

- PE Material strong and durable
- X Ray compatible
- **Add:** Head Immobilizer (DTB-1), Spider strap (DDC-206)

### DIMENSIONS

Length	1830mm
Width	440mm
Height	50mm



## HEAD IMMOBILIZER

DTB-1

### DESCRIPTION

- Material: Polyethylene & PU foam
- Easy to clean and disinfect
- Minimal interference with X Ray / MRI or CT scanning procedures

### DIMENSIONS

N/A



## SPIDER STRAP ADJUSTABLE

DDC-206

### DESCRIPTION

- N/A

### DIMENSIONS

N/A



## SCOOP STRETCHER

DDJ-3B

### DESCRIPTION

- Aluminium construction
- Transfer wounded patients quickly
- Length of stretcher can be adjusted to fit the height of the patient
- Foldable for space saving
- Max weight capacity 159kg

### DIMENSIONS

Extended:

Length	2040mm
Width	440mm
Height	70mm

Folded:

Length	1690mm
Width	440mm



## REQUEST A QUOTE

 [carla@meditek.co.za](mailto:carla@meditek.co.za)

 021 553 1044 | 084 422 6603

 [www.meditek.co.za](http://www.meditek.co.za)