

MEDCOR MEDICAL EQUIPMENT (PTY) LTD t/a



M E D I T E K *hemco*  
MEDICAL EQUIPMENT SALES

# PRODUCT CATALOGUE

## SEATING

---

- > STOOLS
- > CHAIRS
- > RECLINERS
- > SHOWER CHAIR
- > COMMUNE
- > WHEELCHAIR



## ADJUSTABLE HEIGHT STOOL - NO BACKREST

M26

### DESCRIPTION

- Height adjustable stool by means of lever below the seat
- Seat upholstered
- Base chrome finish (without foot ring)
- Mobile on non marking castors

### DIMENSIONS

Diameter	300mm
Height	550 to 850mm



## ADJUSTABLE HEIGHT STOOL - WITH BACKREST

M25

### DESCRIPTION

- Height adjustable stool by means of lever below the seat
- With adjustable backrest by means of knob
- Seat and back upholstered
- Base chrome finish (without foot ring)
- Mobile on non marking castors

### DIMENSIONS

Diameter	300mm
Height	550 to 850mm



## SALLI SWAY - STOOL

SALLIS

### DESCRIPTION

- Divided saddle upholstered seat
- Benefits: Active sitting, Good posture, Healthy back, Good circulation & Deeper breathing
- Hand controlled height adjustment lever
- Rocking mechanism enables the seat to be tilted into every direction without a separate adjustment lever
- 5 Star chrome base
- Chair comes standard with Medium gas cylinder – Max Height 750mm and Lowest 580mm
- **Note: Assembly required**
- **Optional:** Also available in Polyurethane seat

### GAS SPRING SPECIFICATIONS

	Lowest	Highest	User height
Short gas lift	490mm	610mm	Up to 160
Medium gas lift	580mm	750mm	165 – 185
Long gas lift	660mm	915mm	Up to 170



## SALLI CHIN STOOL

SALLIC

### DESCRIPTION

- Divided saddle upholstered seat fixed
- Benefits: Active sitting, Good posture, Healthy back, Good circulation & Deeper breathing
- Hand controlled height adjustment lever
- 5 Star chrome base
- Chair comes standard with Medium gas cylinder – Max Height 750mm and Lowest 580mm
- Max weight capacity 120kg
- **Note: Assembly required**
- **Optional:** Also available in Polyurethane seat

### GAS SPRING SPECIFICATIONS

	Lowest	Highest	User height
Short gas lift	480mm	605mm	Up to 160
Medium gas lift	570mm	740mm	160 – 180
Long gas lift	645mm	910mm	Up to 165



## LBC / GERIATRIC CHAIR

LBCMID - Midback chair  
LBCHI - Highback chair

### DESCRIPTION

- Mild steel construction - epoxy coated finish
- Upholstered vinyl seat and back
- Supawood eating tray (detachable)
- Vinyl colours: Black, Burgundy, Blue or Beige
- Epoxy colours: Black, Brown, White or Beige
- **Optional:** Midback or Highback

### DIMENSIONS

	Midback	Highback
Seat Size	500 x 500mm	500 x 500mm
Seat height (from floor)	470mm	470mm
Back height (from seat)	500mm	750mm



## BO RECLINER

RECBO

### DESCRIPTION

- Upholstered PU leather seat, backrest and arms
- Retractable
- Upholstery in black or brown
- **Note: Assembly required**

### PRODUCT MATERIAL EXPLAINED:

**Leather Touch PU / Polyurethane (PU)** – PU is synthetic leather that is chemically created. PU is known to be a strong, flexible and durable material. It is good for its price and eco friendliness and creates a good leather-like look on furniture. PU furniture is in the budget price range.

### DIMENSIONS

Height	1020mm
Width	760mm
Depth	910mm



## SHOWER CHAIR – HEAVY DUTY

TILTON2

### DESCRIPTION

- Tubular steel construction
- Polypropylene moulded seat
- Seat has 5x 6mm drainage holes
- Height adjustable
- Removable armrests
- Max user weight 165kg
- Mobile with 2x lockable castors
- **Note: Assembly required**



## COMMODOE – HEAVY DUTY

COMDA

### DESCRIPTION

- Height adjustable
- Drop down padded armrests
- Width between armrests 65cm
- Padded soft seat
- Max user weight 450kg
- Non mobile
- Includes a 10L commode bucket
- **Note: Assembly required**



## SHOWER CHAIR / COMMODOE – HEAVY DUTY

ASTON

### DESCRIPTION

- Height adjustable
- Moulded seat for comfort
- Armrest flip up for side transfer
- Adjustable backrest
- Adjustable height footrests / foldaway
- Max user weight 159 kg
- Unit is on lockable wheels
- Includes a 5L commode bucket with lid
- **Note: Assembly required**



## WHEELCHAIR ECONO FOLDING

MEC37

### DESCRIPTION

- Tubular steel construction with black coating
- Fixed arm & footrests
- Black nylon upholstered
- Seat width 45cm
- Foldable
- Max user weight 100kg
- **Note: Assembly required**



## WHEELCHAIR

FREETEC


### DESCRIPTION

- Tension adjustable backrest
- Height adjustable footrest
- Extra strong from castors
- Height adjustable push handles
- Seat width 41cm
- Max user weight 125kg
- Foldable
- **Note: Assembly required**



## REQUEST A QUOTE

 [carla@meditek.co.za](mailto:carla@meditek.co.za)

 021 553 1044 | 084 422 6603

 [www.meditek.co.za](http://www.meditek.co.za)