

MEDCOR MEDICAL EQUIPMENT (PTY) LTD t/a



M E D I T E K *hemco*  
MEDICAL EQUIPMENT SALES

# PRODUCT CATALOGUE

## TABLES & LOCKERS

---

- > OVERBED TABLES
- > BEDSIDE LOCKERS
- > MAYO TABLE



## OVERBED TABLE – ADJUSTABLE HEIGHT

M68

### DESCRIPTION

- Tubular steel construction with baked epoxy coating
- Height adjustable laminated wooden top
- **Note:** Other finishes available

### DIMENSIONS

Top	900x405mm
Height	740-1090mm



## OVERBED TABLE – FIXED HEIGHT

M69

### DESCRIPTION

- Tubular steel construction with baked epoxy coating
- Fixed height laminated wooden top
- **Note:** Other finishes available

### DIMENSIONS

Diameter	910 x 370mm
Height	1020mm



## OVERBED TABLE

MEC33

### DESCRIPTION

- Tubular steel frame with epoxy coating
- Top sheet metal with epoxy finish
- Fixed height
- **Note** Wooden top & S/steel top also available

### DIMENSIONS

Length	1030mm (to fit over width of bed)
Width	380mm
Height	1015mm



## BEDSIDE LOCKER MELAMINE

M91

### DESCRIPTION

- Melamine finish – Oak
- One telephone drawer
- One standard drawer
- **Note :** More finishes available
- **Optional:** Two doors, one with a shelf and one with a rail

### DIMENSIONS

Height	950mm
Width	550mm
Depth	575mm



## BEDSIDE LOCKER MELAMINE

M92

### DESCRIPTION

- Melamine finish – Oak
- One door with internal shelf
- Two drawers - One of which is a telephone drawer
- Magazine rack
- **Note:** More finishes available

### DIMENSIONS

Height	800mm
Width	490mm
Depth	490mm



## BEDSIDE LOCKER

M72

### DESCRIPTION

- Sheet metal construction – mild steel with baked epoxy coating
- S/steel top
- Upper section with one drawer
- Lower section with one door and internal shelf
- Fitted with a towel rail
- **Optional:** With wooden top

### DIMENSIONS

Height	850mm
Width	430mm
Depth	430mm



## BEDSIDE LOCKER

MEC12 - Complete epoxy coated  
MEC12A - S/steel top

### DESCRIPTION

- Sheet metal construction with baked epoxy coating
- Epoxy coated top and body (MEC12)
- Epoxy coated body with s/steel top (MEC12A)
- Upper section with one drawer (pull-out)
- Lower section with one door and internal shelf

### DIMENSIONS

Height	850mm
Width	430mm
Depth	430mm



## BEDSIDE LOCKER - COMBO

A2-032-EF2

### DESCRIPTION

- Mild steel construction with epoxy coated finish
- S/steel tops
- Bedside locker with integrated 950mm H wardrobe section
- Locker section with cabinet and one drawer
- With towel rail

### DIMENSIONS

Length	720mm
Width	430mm
Height	1080mm

### CASTORS:

6 x 50mm swivel castors



## MAYO TABLE

MEC5EPY  
MEC5SS

### DESCRIPTION

- MEC5EPY - Tubular mild steel epoxy coated frame
- MEC5SS - Tubular s/steel frame
- With removable stainless steel tray - 600mm x 400mm
- Height adjustable 960mm - 1270mm

### DIMENSIONS

Length	600mm
Width	400mm
Height	960 - 1270mm



## MAYO TABLE

M58 / M58SS

### DESCRIPTION

- Tubular steel frame – epoxy coated
- Stainless Steel tray
- Height adjustable – 960mm – 1270mm
- **Optional:** Complete s/steel unit (M58SS)


### DIMENSIONS


Length	600mm
Width	400mm
Height	960-1270mm



## REQUEST A QUOTE

 [carla@meditek.co.za](mailto:carla@meditek.co.za)

 021 553 1044 | 084 422 6603

 [www.meditek.co.za](http://www.meditek.co.za)