

MEDCOR MEDICAL EQUIPMENT (PTY) LTD t/a



M E D I T E K *hemco*  
MEDICAL EQUIPMENT SALES

# PRODUCT CATALOGUE

# TREATMENT EXAMINATION

# TABLES

- > COUCH
- > WOODEN COUCH
- > UNIVERSAL COUCH
- > ELECTRIC COUCH
- > NEURO PLINTH



## EXAMINATION COUCH

MEC20

### DESCRIPTION

- Mild steel tubular construction with epoxy coated finish
- Fully upholstered top
- With ratchet adjustable backrest
- **Optional:** Breathing slot
- **Optional:** Upholstery can be done in White, Navy blue or Black
- **Optional:** Epoxy can be done in White or Ivory

### DIMENSIONS

Length	1830mm
Width	610mm
Height	760mm



## EXAMINATION COUCH - WOODEN

M6A

### DESCRIPTION

- Wooden construction (Meranti)
- Fully upholstered top
- With ratchet adjustable backrest

### DIMENSIONS

Length	1830mm
Width	610mm
Height	760mm



## COUCH - UNIVERSAL EXAMINATION

M96

### DESCRIPTION

- Tubular steel construction with epoxy finish (also available with Trendelenburg or Stainless Steel)
- Adjustable leg and head sections
- Upholstered with vinyl
- With holder for stainless steel bowl
- With two Sheppard hooks
- **Optional:** Stainless Steel Bowl

### DIMENSIONS

Length	1800mm
Width	560mm
Height	780mm



## 2 SECTION COUCH – ELECTRIC EXAMINATION

2SECELEC

### DESCRIPTION

- 2 section table
- Gaslift backrest with positive / negative tilt
- Electric height by means of foot switch
- With breathing slot
- Retractable wheelbase
- Max weight capacity 240kg

### SPECIFICATION

Head Section Min Height	40° Negative
Head Section Max Height	90° Positive
Minimum Height	450mm
Maximum Height	9200mm
Total Width	700mm
Total Length	2000mm



## 3 SECTION COUCH – ELECTRIC EXAMINATION

3SECEL

### DESCRIPTION

- 3 section table
- Gaslift backrest with positive / negative tilt
- Electric height by means of foot switch
- Adjustable leg section by means of gaslift
- With breathing slot
- Retractable wheelbase
- Max weight capacity 240kg

### SPECIFICATION

Minimum Height	450mm
Maximum Height	920mm
Head Section Max Height	90°
Head Section Min Height	40°
Leg Section Max Height	30°
Leg Section Min Tilt	25°
Total length	2000mm
Total Width	700 mm



## NEURO PLINTH 2 SECTION

NEURPLIN2


### DESCRIPTION

- Electric height adjustable
- 6 cm thick foam padding
- Durable vinyl upholstery
- Mild steel epoxy powder coated for easy cleaning and disinfecting
- Fully retractable wheelbase for ease of movement
- Head section adjustment with two gas struts allowing for one hand operation
- Electric height adjustable by means of footswitch
- Lifting Capacity of 330Kg
- Breathing slot with plug



## REQUEST A QUOTE

 [carla@meditek.co.za](mailto:carla@meditek.co.za)

 021 553 1044 | 084 422 6603

 [www.meditek.co.za](http://www.meditek.co.za)